




COVID-19 Antibody Tests (CA)

Product Specification Sheet



INSTI® COVID-19 Antibody Tests (50 individual tests with support materials)

Product Code	90-1092
List Price	\$
Sample Type	Human venipuncture whole blood, serum and plasma

Kit Contents	Materials: 50 individual tests per box
	1 Single Test with Support Materials <ul style="list-style-type: none">• 1 Membrane Unit• 1 Bottle 1, Sample Diluent (1.5 mL)• 1 Bottle 2, Colour Developer (1.5 mL)• 1 Bottle 3, Clarifying Solution (1.5 mL)• 1 Single-Use Safety Sterile Lancet• 1 Capillary Test Pipette• 1 Alcohol Swab• 1 Instructions for Use

Test Run Time	1 minute
---------------	----------

Shelf Life ¹	6 months
-------------------------	----------

QC Requirements	<ul style="list-style-type: none">• For new INSTI® operator verification prior to performing testing on patient specimens• When switching to a new lot number of INSTI test kits• Whenever a new shipment of kits is received• When temperature during storage of the kit falls outside of 2 - 30°C (35.6 - 86°F)• When the temperature of the test area falls outside of 15 - 30°C (59 - 86°F)• Approved for lab use only
-----------------	---

Storage Conditions	2 - 30°C (35.6 - 86°F) until expiration date
--------------------	--

Performance ²	Sample Type	PPA	NPA
	Plasma/Serum	98.2% (163/166)	99.6% (259/260)

¹From date of manufacture. Shelf life will be updated/extended as further stability data becomes available.

Performance estimated using SARS-CoV-2 PCR positive plasma/serum samples and negative plasma/serum samples, with negative samples collected October 2019 or prior.

Additional Accessories (Sold Separately)	Code	Product
	10-1069	Lancets, 200 pack
	80-1087	Test Capillary Pipettes, 48 pack

Manufactured by:



bioLytical Laboratories Inc.

406-13251 Delf Place
Richmond, BC Canada V6V 2A2
1-866-674-6784 | info@biolytical.com

insti.com



65-4017(A) 01/2022